

Measuring Temperature

Aim: At the end of the training session participants Should perform the steps in order and correctly

Tools Thermoscan

Steps:

- 1.Wash your Hands
- 2.Explain the procedure to the patient
- 3.Handle the thermoscan and control the screen.(around 34 C)
4. Insert the thermoscan to the pprobe (disposable)
5. Place the probe in the ear but do not push too hard.
- 6.Press the button which runs the thermoscan, wait until hearing the electronic alert
- 7.Retrack the thermoscan and read the value on screen
- 8.Record the result, inform the patient and wash your hands.